

SET UP ON SMALL OPTI - BENCH. USE AERO EKTAR LENS AND SMALL LIGHT SOURCE.

PLACE A SLIDE IN THE SLIDE HOLDER

PLACE CONDENSERS BACK TO BACK. PLACE CONDENSERS AND SLIDE HOLDER AS CLOSE TOGETHER AS POSSIBLE.

PLACE CONDENSERS AND SLIDE HOLDER 145 mm FROM THE LIGHT SOURCE.

PLACE AERO EKTAR LENS 420 mm FROM THE LIGHT SURCE.

FOCUS IMAGE ONTO SCREEN BY MOVING CART CLOSER OR FARTHER AWAY.